

Silver Dressings

Silverlon® IV/Catheter Dressings are sterile antimicrobial barrier dressings designed to reduce bacterial colonization associated with catheter related infections.

Applicable to: Peripheral IV Catheters, Central Venous Lines, Arterial Catheters, Peritoneal Dialysis Catheters, Dialysis Catheters, External Fixator Pins, and Drain Tubes.



Silverlon® Tubular Stretch Knit Products are designed to bring antimicrobial protection to digits and certain hard to dress extremities. These easy to use products are gently stretched to apply, and can be used as primary dressings. The Silverlon® dressings should then be covered or secured with conventional technique.

The **Silverlon® Burn Wrap** is a broad spectrum 7 day use antimicrobial barrier dressing. The surface of the nylon-based dressing is plated with pure metallic silver, which is the source of the antimicrobial silver ion. The silver ions disassociate from the surface of the dressing in the presence of moisture. Silver ion release rates are maintained by maintaining the moisture content of the burn and wound dressing.

Silverlon® Acute Burn Gloves are open finger and oversized 20% for use immediately after acute injury or post surgery for both partial to full thickness burns. They offer broad antimicrobial protection to often hard to dress parts of the hand.



Silverlon® Wound Contact Dressing is safe and effective for any laceration, abrasion, cut, burn, or chronic wound. Even works on most rashes.

Effective Treatment for stubborn MRSA Infections. Delivers a speedy and sustained kill against MRSA Bacteria.

When kept clean and moist it provides powerful antimicrobial and anti inflammatory action for up to 7 days.



Silverlon® Surgical Wound Dressing (Island Dressing) can be used effectively to treat, and heal virtually any primary incision site in the Operating Room, and in the Hospital prior to patient discharge. The Silverlon® Surgical Wound Dressing offers advanced antimicrobial protection for the following surgical procedures: Cardiac Surgery (Cardio Vascular Surgery), Colon Surgery, Abdominal Hysterectomy, Abdominal Caesarean (C Section), Vascular Surgery - Repair of Aneurysm, Vascular Surgery - Thromboendarterectomy, Vascular Surgery - Vein Bypass, Hip Arthroplasty, and Knee Arthroplasty Surgery. The Silverlon® Surgical Dressing is applicable for treatment of surgical incisions, partial thickness burns, skin grafts, donor sites, lacerations, abrasions, dermal ulcers (vascular, venous, pressure and diabetic).



Silverlon® Negative Pressure Wound Dressings are applied for improved outcomes with Negative Pressure Wound Therapy Systems when treating chronic and acute wounds, burns, and surgical procedures.

The Silverlon® NPD Advantage:

- 7 day use over skin grafts while under negative pressure therapy systems (1 dressing kit for 1 week of therapy)
- Fewer dressing changes over the course of therapy
- Produces an accelerated wound healing potential
- No infections reported
- Provides medical and economical advantages
- Minimal reduction in fluid flow capacity as compared to VAC foam alone

