

# AME 1400 LOW AIR LOSS MATTRESS



## AME1400

The AME1400 is a low air loss mattress overlay system that revolutionizes the overlay market by providing customized comfort and optimal skin protection! The AME1400 low air loss systems with the ultimate skin protection of a ROHO® HIGH PROFILE® Cushion in the heel area to create an optimal environment for skin healing and protection.

The AME1400 is so versatile! It sets up in minutes on top of any existing mattress – no need to store mattresses! The ROHO HIGH PROFILE cushion in the heel area, our HealFloat™ feature, protects heels and creates an optimal environment for wound healing. Simply remove the cushion to use on any chair for added skin protection.



- System Lock: Helps to prevent accidental adjustments, and automatically resets the system in case of power failure.
- Firm Mattress Mode: Helps facilitate nursing care during dressing changes and patient repositioning.
- HealFloat™: ROHO HIGH PROFILE Cushion inserted in the heel area creates an optimal environment for skin healing and protection.
- Pulsation (alternating pressure): Provides two cycle settings that readjust the system to ensure patient comfort.



### Key Product Features

- Fits on top of existing mattress. No need to store mattress!
- Easy-to-adjust control unit provides constant inflation.
- Adjustable control settings provide optimum comfort.
- Mattress overlay and cushion are easy to clean.
- Durable, high-grade materials provide long product life.
- Low air loss overlay system
- Air cells are organized in three zones (heat, seat, foot)
- ROHO high profile cushion is integrated at the foot section
- Integrated side air bolsters to prevent entrapment and keeps patient centered
- Low Air loss circulates air for moisture management
- Waterproof, breathable high quality nylon top cover reduces friction and shear
- Optional foam supplement can be used in place of existing mattress

# AME 1400 LOW AIR LOSS MATTRESS



## AME1400

### Easy-to-operate

This powerful low-air-loss system gives you advanced, easy-to-operate settings in a simple overlay. With the help of AME1400, it's never been this easy to protect patients' skin and keep them comfortable.

### Powerful

Control unit provides constant inflation and can be adjusted to suit your needs.

### Zoned Overlay

The AME1400 low air loss mattress overlay system creates an optimal environment for skin healing and protection. The zoned overlay utilizes an internal manifold which allows three independent mattress sections to be set using one button on the control unit. AME1400 manages moisture and temperature in order to promote skin healing with "true" low air loss allowing air to escape through the top of the mattress. Comfortable Adjustable control unit settings for optimum comfort.

### Durable

High grade materials provide long product life even if patient's needs change.

### Firm Mattress Mode

To aid in the placement and maneuvering of the patient, the control unit is equipped with a Firm Mattress Mode. The Firm Mattress Mode will fully inflate the mattress making it firm for easy positioning. The Firm Mattress Mode will automatically terminate when the time displayed on the keypad has expired.

### Pulsate Mode

This control unit is equipped with the Pulsate Mode. When selected, the control unit will continuously readjust the system pressure (proportional to the Comfort Level setting). Change in pressure is gradual to ensure maximum patient comfort.

