

# AME 3500 AND AME 3584

## LOW AIR LOSS AND ALTERNATING PRESSURE MATTRESS SYSTEM.



### AME3500 / AME3584

The American Medical Equipment AME3500 and AME3584 system provides alternating pressure and is designed to be used in the prevention, treatment and management of Pressure Ulcers. Pump produces consistent air flow. Special vibration isolating features insure that the patient rests easily, without feeling the pumps vibrating action. LED and audible alarm provide critical performance data and alert caregiver when pressure falls below 10mmHg (inner mattress pressure). Easy to replace individual bladders. Mattress has 20 air bladders, nine (9) with laser holes to achieve low air loss function, to increase patient comfort, keep the patient cool and speed recovery. Spring lined air tubes prevent kinks and disruption of air flow. CPR valves provide fast (20-Second) deflation.

- Quiet, portable, lightweight and very reliable, the SelectAir mattress and control unit provide a cost-effective alternative to other pressure management systems.
- Firm Mattress Mode: Helps facilitate nursing care during dressing changes and patient repositioning.
- Upright Patient Mode: Adds support when the head of the bed is elevated.
- Foot Pillow: Helps to prevent unnecessary skin breakdown due to patient migration.



### Key Product Features

- Air cells are organized in three zones (head, seat, foot) .
- Integrated side air bolsters to prevent entrapment and keeps patient centered.
- Low air loss circulates air for moisture management.
- Alternating pressure provides three cycle settings continually readjusting the system.
- Auto firm to assist in transfers and repositioning.
- Fowler (upright patient) supports patient when head of the bed is elevated.
- Auto comfort mode simplifies setup and is based on patient's weight distribution.
- Waterproof, breathable high quality nylon top cover reduces friction and shear.
- Available with integrated raised perimeter with egress openings for fall prevention.
- Available in 84" Length (AME 3584) for taller patients
- Low air loss and alternating pressure mattress system.

# AME 3500 AND AME 3584

## LOW AIR LOSS AND ALTERNATING PRESSURE MATTRESS SYSTEM.



### AME3500 / AME3584

#### User-Friendly Control Unit

Rapid inflation, easy set-up, and a motor so quiet, the patient will hardly know the Quick-Change Filter is there!

#### Integrated Side Bolsters

Keeps patient centered in the bed and reduces the risk of entrapment.

#### Keeps Patient Cool and Dry

"True" Low Air Loss technology allows air to escape from the top of the mattress and directly interact with moisture at the patient/ bed boundary.

#### Wicking Action

Air escapes through tiny holes on the top of the mattress, and circulates under a vapor-permeable cover, to wick moisture away from patient.

#### Exceptional Foot Protection

Foot pillow to address patient migration.

#### Quick-Change Filter

Filter is easy to access, and quick to change for hassle-free maintenance.

