

Dome Pad with double lumen design feature and benefits:

1. It will prevent blockage alarming, if there is a blockage at one of the tubing, the other tubing will deliver negative pressure to the wound site. Less trouble-shooting and less user frustration.
2. The double lumen was designed to prevent exudate pooling and wound maceration, there is constant negative pressure in the wound site with the double lumen design.
3. The ridge ceiling design inside the dome, is guaranteed not to collapse even if the dome is placed in an area where there is tremendous pressure.
4. The dome pad skirt drape was created to provide easy placement and proper sealing.
5. The dome pad connector is universal and can be attached to any OEM npwt pump system.
6. The dome and tubing will disconnect, in an event of a pulling or equipment dropping. The dome will remain attached to the drape and prevent infection control issue.
7. The tubing is very soft to the skin and will not create any shearing or friction to the patient skin.

