
ALTERNATING PRESSURE OVERLAY SYSTEM AME300



AME 300 MATTRESS OVERLAY

Designed using deep cell air technology for international use, the AME300 Mattress has 20 nylon-coated, 8-inchhigh air cells to immerse and envelop the patient for optimal pressure redistribution. The AME300 Mattress also features a foam base for additional protection during power outages and is covered with a quilted Plexus[®] nylon top cover that reduces shear. Complete with an easy to use AME300 Controller featuring a pressure selection dial control, the AME300 Mattress System is an ideal pressure management solution for home care and healthcare facilities.



Key Product Features

- Selectable therapy mode of reactive or alternating pressure
- Alternating pressure every five minutes for enhanced therapy
- Modular design facilitates cell interchangeability, cleaning and set-up
- Provides pressure management throughout entire alternating cycle
- Plexus[®] quilted top sheet is highly vapor permeable and oversized, providing low friction and low shear. Machine washable, the cover is available in Nylon or Silver3[®], which is embedded with silver ions for added protection
- Quick hose release for rapid CPR deflation
- Compatible with the optional Foot Stop/ZIP Cushion[™] to help prevent foot drop while floating the heel to eliminate heel pressure
- Lightweight and compact for ease of transporting
- User-friendly control unit allow the caregiver to adjust for customized therapy;
- Equipped with "Max. Inflate" for easy patient positioning;
- Patient settings are restored after power outages for nursing convenience;
- Lightweight and compact for easy transport; and
- CPR quick-release for rapid deflation.

Cycle time: 5 minutes. Weight capacity: 350 lbs. Dimensions: 34"w x 76"l x 5"h. Warranty: 2-yr limited on air controller and 1-yr on air mattress.



Phone: 1-866-633-2337 Fax: 1-614-237-1177 customerservice@ame-medical.com