



## AME PRESSUREGUARD® APM<sup>2</sup>

The **PressureGuard® APM<sup>2</sup>** enables you to change instantly and easily from alternating pressure to basic lateral rotation with a single, recessed toggle switch. The silent, reliable, maintenance-free control unit allows simple adjustment for user comfort.

With its two distinct therapy modes, the APM2 gives you twice the treatment options. Twice the flexibility. Twice the value. Change instantly and easily from alternating pressure to basic lateral rotation with a single, recessed toggle switch. Requires no disconnecting or exchanging of air lines. And, when ordered with the new Deluxe control unit, the APM2 offers even more value. Selectable audible alarm can be engaged to alert caregiver to low pressure status. The "Low Pressure" indicator light gives users at-a-glance assurance of proper function.



### Key Product Features

- Therapeutic Geo Matt anti-shearing top surface
- Two-part Safety Edge; Smooth, subtle shifts in weight-bearing in ten-minute cycles without sharp pressure drop-offs
- Weight limit 350 lbs.
- Mattress weight 20 lbs.
- Pump weight 5.25 lbs.
- 18 Month warranty. UL listed

#### Primary Use:

- Aggressive Treatment of Pressure Ulcers

#### Patient Indications:

- Cant or wont reposition frequently
- Can be positioned off of ulcer in at least 2 positions
- Prevention for those at high risk for skin breakdown(e.g. Braden score 12 or less)

#### Key Advantages:

- Two therapies in one mattress
- Single switch operation
- PressureGuard® design: stability, safety, anti-shearing

